

SPELEC

The intuitive spectroelectrochemistry

PEOPLE YOU CAN TRUST



The best of two worlds

Spectroelectrochemistry is a hyphenated technique that takes the advantages of electrochemistry and optical methods in a single experiment, providing the whole vision of chemical processes that take place on the electrode surface. For this reason, spectroelectrochemistry is currently one of the most popular techniques for gathering molecular, kinetic, and thermodynamic information from the reactants, intermediates, and/or products involved in electron transfer processes being useful in a huge variety of fields.

ONE SOFTWARE FOR

ADVANCED RESEARCH

In comparison with the modular standard

run spectroelectrochemical measurements.

DropView SPELEC is the only dedicated spec-

data treatment tools making accessible time-

resolved spectroelectrochemistry to everyone.

troelectrochemical software with specific

set-ups, only one software is required to

SPELEC instruments, based on the operando concept, are the only fully- integrated equipments dedicated to spectroelectrochemistry. All components, (bi)potentiostat/galvanostat as well as light source and spectrometer are integrated in the same instrument. This configuration allows that a key factor, such as real synchronization (non-triggered), is achieved between both techniques.

«Spectroelectrochemistry is a powerful technique that combines electrochemistry and spectroscopy»



Spectrometer connector

Available models

	SPELEC	SPELEC1050	SPELECNIR	SPELECRAMAN	SPELECRAMAN638	SPELECRAMAN532
LIGHT SOURCE	200-400 nm (deuterium) 400-1100 nm (halogen)	200-400 nm (deuterium) 400-1100 nm (halogen)	400-2500 nm (tungsten halogen)	785 nm laser	638 nm laser	532 nm laser
SPECTROMETER	200-900 nm	350-1050 nm	900-2200 nm	50-3000 cm ⁻¹	60-4350 cm ⁻¹	70-4500 cm ⁻¹
POTENTIOSTAT	\pm 4 V \pm 40 mA	\pm 4 V \pm 40 mA	\pm 4 V \pm 40 mA	\pm 4 V \pm 40 mA	\pm 4 V \pm 40 mA	\pm 4 V \pm 40 mA
DIMENSIONS	25 x 24 x 11 cm (L x W x H)	25 x 24 x 11 cm (L x W x H)	25 x 24 x 11 cm (L x W x H)	25 x 24 x 11 cm (L x W x H)	25 x 24 x 11 cm (L x W x H)	25 x 24 x 11 cm (L x W x H)

3 INSTRUMENTS IN ONLY 1

Although SPELEC instruments are designed for performing spectroelectrochemical measurements, they can be also used independently and activated/deactivated via software only as (bi)potentiostat/ galvanostat, or as an optical instrument.

COVERING ALL RANGES

SPELEC instruments cover a wide spectral range for performing UV-VIS, VIS-NIR, NIR or Raman spectroelectrochemistry. Depending on your system under study, and the ranges of interest, you may choose between the different compact fully integrated models available.

Just one software dedicated for spectroelectrochemistry...

SPELEC instruments are controlled by DropView SPELEC, the only software in the market dedicated to spectroelectrochemistry. It is thought and designed for (and by) those that want to treat spectroelectrochemical data in just one click.

ANALISYS IN ONE CLICK

Specific spectroelectrochemical data treatment tools are available to facilitate the analysis. Discover the one-click functions accessible to new and expert users.



«Real synchronization matters»

REAL SYNCHRONIZED DATA

Electrochemical and spectroscopic signals are independent (non-triggered) but simultaneous. Optical acquisition starts/finishes measuring at the same time than the electrochemical reaction, obtaining information during the whole experiment.

LIVE ACQUISITION

Spectra are continuously recorded and shown during the measurement. You do not need to wait until the end of the experiment to understand the behavior of the system under study.

AUTOMATIC CONTROL

The shutter of the lamps is also automatically controlled for setting dark and reference spectra as well as for running the experiment. In addition, the laser power is controlled with DropView SPELEC, allowing you optimize it according to the properties of the sample under study.

... with many data treatment tools.



EXPERIMENT FILM

Overview your experiments and save them as video files.



VISIBLE SPECTRUM AT EC POINT

Easy visualization of the spectrum associated with each electrochemical point.



3D PLOT

Outstanding representation of the data and easy visualization of the results.



SPECTRA VS EC

Track the evolution of the optical signal with potential.



DATATABLE

Export straightforward EC and optical data, showing their perfect synchronization.

WORKSPACE

3 axis representation to compare optical and electrochemical signals in a single graph.



BASELINE CORRECTION

Remove the background signal to define the relevant bands.

	EC					Ор	tical		
t(s)	E (r)	1 (uA)	1(8)	36.28L	414.91L	765.61L	1091.22L	1394.20L	1676.68L
0.000	0.000000	-3.193750	0.000	15.820816	693.1290	730.4402	821.8574	773.4597	537.5136
0.100	0.002000	-3.043750	0.100	15.820816	693.1290	730.4402	821.8574	773.4597_	537,5136
0.200	0.004000	-2.907917	0.200	15.820816	693.1290	730.4402	821.8574	773.4597_	537.5136
0.300	0.006000	-2.785833	0.300	15.820816	693.1290_	730.4402	821.8574	773.4597_	537.5136
0.400	0.008000	-2.677083	0.400	15.820816	693.1290_	730.4402	821.8574	773.4597	537.5136
0.500	0.010000	-2.571667	0.500	15.820816	693.1290	730.4402	821.8574	773.4597	537.5136
0.600	0.012000	-2.472500	0.600	15.572894	693.6877_	731.4695	819.3914	773,2119	534,9128
0.700	0.014000	-2.376667	0.700	15.324971	694 2465	732.4988	816.9253	772 9540	532 3121
0.800	0.016000	-2.285417	0.800	15.077049	694 8053	733.5281	814.4592	772,7161	529.7114
0.900	0.018000	-2.198333	0.900	14.829127	695 3640	734 5575	811,9932	772.4683	527.1107
1,000	0.020000	-2.117083	1,000	14.581204	695 9228	735 5868	809 5271	772 2204	524 5100
1.100	0.022000	-2 037500	1.100	14703064	695.5741	735.9775	8097834	771.8045	524 8335
1,200	0.024000	-1962500	1200	14 824925	695 2254	734 3482	810 0396	771 3885	626 1570
1,300	0.026000	-1885417	1,300	14 946785	694 8767	736,7589	810,2959	770 9726	525 4805
1.400	0.028000	-1.015033	1.400	15.069645	694 5290	737.1496	810 5522	770.5567	525.0040
1.500	0.030000	-1747500	1.500	15.190505	694 1792	737 5403	810.8084	770 1408	626 1275
1.600	0.032000	1684583	1600	15 393345	694 8523	739 1545	8110781	772 6950	527 2712
1700	0.034000	-1622500	1,700	15.595585	695 5253	740 7686	811 3478	775 2511	528.4148
1,800	0.036000	-1.562917	1,800	15,798124	696.1984	742 3827	811.6176	777.8063	529 5585
1,900	0.038000	-1505417	1,900	15.000554	695 8714	743 9968	811.8873	780 3614	530 7021
2,000	0.040000	-1.445557	2,000	15 203204	697 5445	745.6109	\$12 1570	782 9166	531.8458
2 100	0.042000	-1 390833	2 100	16 349436	7010517	744 1439	812 7065	780 2396	533 0003
2,200	0.044000	-1340417	2,200	15.495668	704 5589	742 6768	813 2560	777 5626	534 1549
2 300	0.046000	-1292500	2 300	15.641899	708 0661	741 2097	813 8055	774 8855	535 3095
2,400	0.048000	-1247083	2,400	15,788131	7115733	739 7426	814 3550	772 2086	535.4540
2 500	0.050000	-1 109167	2 500	15 934363	715.0805	738 2765	814 9045	769 5316	637 6186
2 600	0.052000	1154583	2 600	15 228415	7110137	739 1830	\$16 2153	772 5901	537 2489
2,700	0.054000	-1.108333	2,700	15.522468	706 9469	740.0905	817.5261	775.6485	536.8791
2,800	0.056000	-1.053750	2.800	14,816521	702 8801	740 9980	818 8368	778 7070	536 5094
2,900	0.058000	-1021667	2 900	14,110573	698 8133	741 9054	820 1476	781 7655	536 1397
3,000	0.050000	0 980417	3,000	13 404526	6947454	742 8129	821 4583	784 8239	535 7699

Looking for spectroelectrochemical accessories?

SPELEC instruments can be used with any kind of cells and set-ups, however you have also available Metrohm DropSens cells and accessories specifically designed for the SPELEC line.

Simplify your spectroelectrochemical set-up with optical cells to perform reflection or transmission experiments with SPEs or conventional electrodes. Complete your optical configuration with the fibers, probes and accessories from our catalogue.



SPECTROELECTROCHEMICAL CELLS FOR SCREEN-PRINTED ELECTRODES

Cells suitable to perform UV-VIS, NIR or Raman spectroelectrochemical experiments with screen-printed electrodes. Reflection and transmission cells are available. Benefit from an easy sensor replacement thanks to an opening and closing system with magnets.



SPECTROELECTROCHEMICAL CELLS FOR CONVENTIONAL ELECTRODES

Cells suitable to perform UV-VIS, NIR or Raman spectroelectrochemical experiments with conventional electrodes. Reflection (UV-VIS, NIR and Raman) and transmission (UV-VIS and NIR) cells are available, allowing an easy optimization of the focal distance. If your spectroelectrochemical setup includes a microscope, a specific Raman cell is designed for these measurements.



SPECTROELECTROCHEMICAL REFLECTION CELLS FOR THIN-LAYER FLOW-CELL SCREEN-PRINTED ELEC-TRODES

Cells for performing UV-VIS, NIR or Raman spectroelectrochemical measurements in flow conditions in combination with the TLFCL-CIR format SPEs.



REFLECTION AND TRANSMISSION UV-VIS AND NIR FIBERS

Optical fibers designed for UV-VIS or NIR measurements. Reflection probe shows 6-around-1 fiber bundle design, allowing the optimization of the optical measurements.

Ordering codes and description

UV-VIS and NIR	
DRP-CLENS	Collimator lens for TRANSCELL
DRP-CUV	Cuvette holder for PTGRID-TRANSCELL
DRP-FLKIT	Fluorescence kit
DRP-FLKITSPE	Fluorescence kit for screen-printed electrodes
DRP-LEDRGB	LED light red green blue
DRP-LEDUV275	LED light-UV 275 nm
DRP-LEDVIS395	LED light-VIS 395 nm
DRP-PTGRID-TRANSCELL	Transmission cell with conventional electrodes
DRP-REFLECELL	Reflection cell for screen-printed electrodes
DRP-REFLECELL-C	Reflection cell for conventional electrodes
DRP-REFLEPACK-VIS-UV	Pack for reflection experiments with screen-printed electrodes
DRP-RPROBE-VIS-NIR	Reflection probe VIS-NIR
DRP-RPROBE-VIS-UV	Reflection probe VIS-UV

DRP-TFIBER-VIS-NIR	Transmission fiber VIS-NIR
DRP-TFIBER-VIS-UV	Transmission fiber VIS-UV
DRP-TLFCL-REFLECELL	Reflection cell for thin-layer f
DRP-TRANSCELL	Transmission cell for transpar
DRP-TRANSPACK-VIS-UV	Pack for transmission experir
DRP-VKITSPELEC	Verification kit for SPELEC an
DRP-VKITSPELECNIR	Verification kit for SPELECNIF
RAMAN	
DRP-RAMANCELL	Raman cell for screen-printed
DRP-RAMANCELL-C	Raman cell for conventional
DRP-RAMANCELL-M	Raman cell for microscope
DRP-RAMANPROBE	Raman probe for SPELECRAN
DRP-RAMANPROBE532	Raman probe for SPELECRAN
DRP-RAMANPROBE638	Raman probe for SPELECRAN
DRP-TLFCL-RAMANCELL	Raman flow-cell for thin-laye
DRP-VKITSPELECRAMAN	Verification kit for SPELECRA
DRP-VKITSPELECRAMAN532	Verification kit for SPELECRA
DRP-VKITSPELECRAMAN638	Verification kit for SPELECRA



RAMAN PROBES FOR DIFFERENT WAVELENGTHS

Raman probes allow the excitation of the sample and the collection of Raman signal. Different wavelength probes are available: 532, 638 and 785 nm.

flow-cell integrated screen-printed electrodes

rent screen-printed electrodes

ments with screen-printed electrodes

nd SPELEC1050

2

d electrodes

electrodes

MAN

MAN532

MAN638

er flow-cell integrated screen-printed electrodes

MAN

MAN532

MAN638

Why spectroelectrochemistry will improve your research?

Spectroelectrochemistry provides valuable information in a huge variety of different fields. Discover more applications!





ELECTROCATALYSIS

- Analysis of the structure and the electrocatalytic activity of new catalysts
- Monitoring and quantification of water oxidation reaction, oxygen and hydrogen evolution, reduction of CO₂, ammonia oxidation, etc.
- Determination of the stability and conversion of complexes



MATERIAL SCIENCE

- Characterization of morphological and structural properties of carbon nanostructures, nanoparticles, polymers, perovskites, semiconductors, nanocrystals, alloys, composites, etc.
- Study and monitoring of electronic states and doping and degradation processes
- Evaluation of electrochromic capabilities and modulation of optical properties



FUNDAMENTAL CHEMISTRY

- Elucidation of reaction mechanisms
- Identification and quantification of intermediates and generated products
- Calculation of electrochemical and optical parameters

ORGANIC AND INORGANIC CHEMISTRY

- Investigation of vibrational properties of phthalocyanines, porphyrins, organometallics, coordination complexes, dyes, etc.
- Study of the stability, degradation and redox-induced color changes
- Understanding electronic coupling, charge delocalization, intraligand transitions and oxidation ligand-to-ligand charge transfer transitions



CORROSION

- Identification and monitoring of the transformation of products generated during corrosion processes
- Study of the inhibition of corrosion processes on metallic (copper, iron, etc.) surfaces
- Elucidation of the interactions and the reactions of different ions with metallic surfaces





ENERGY

- Elucidation of the structure-property relationships of material employed in solar cells.
- Study of the degree of charge delocalization in rechargeable batteries



SENSING

- Detection and quantification of a huge variety of chemical compounds
- Development of novel sensing platforms to improve the sensitivity capabilities
- Overcoming the limitations of other methods, avoiding complex instrumentation and tedious protocols, saving cost, time, etc.

BIOLOGY AND LIFE SCIENCE

- Monitoring of denaturation, renaturation, hybridization, and interaction processes
- Detection of biological compounds and resolution of biological mixtures
- Study of DNA, neurotransmitters, antioxidants, antitumor agents, proteins, enzymes, etc.





ENVIRONMENTAL

- Monitoring of pollutants in filtration processes
- Quantification of painkillers in wastewater
- Direct observation of amalgamations

Quality



Further information

Metrohm DropSens is a company certified in ISO 9001 and in ISO 13485 (for the 'manufacture of sensors for medical devices') Quality Management Systems Please contact your Metrohm representative or Metrohm DropSens at sales.dropsens@metrohm.com

Metrohm DropSens, Parque Tecnológico de Asturias C/ Faya 28 33428 Llanera (Asturias) Spain www.metrohm-dropsens.com

www.metrohm.com