



Insplorion

# Sensors



- *Real-time information about refractive index changes at surfaces.*
- *Study the influence of surface chemistry.*
- *Asses if and how controlled topography/surface structure influences molecular adsorption.*

## Insplorion™ sensors

Insplorion provides NanoPlasmonic Sensing (NPS) sensors, which enable highly reproducible ultrasensitive measurements of refractive index changes close to (< 30 nm from) the sensor surface. The sensors can be coated with a range of materials, which allows studies of how surface chemistry influences processes such as molecular adsorption and thin film phase transitions. Insplorion also offers plasmonic nanostructures of various types and sizes, allowing the user to perform systematic studies of how controlled topography/surface structure may influence surface processes.

**Insplorion™**

Sahlgrenska Science Park


Medicinaregatan 8A

413 90 Göteborg

Sweden

[www.insplorion.com](http://www.insplorion.com)

## XNano/X1 sensors




The dimensions of XNano and X1 sensors are 9.5x9.5x1 mm. The NPS sensing structure is deposited on a fused silica substrate. The sensor surface can be covered by a top coating if desired.

The sensing structures are distributed evenly over the surface.



## Acoulyte sensor



Acoulyte sensors are QCM-D sensors with a SiO<sub>2</sub> spacer layer deposited on the top gold electrode. The NPS structure is then deposited on top of this layer. The sensor surface can be covered by a top coating if desired.

The diameter of the Acoulyte sensors is 14 mm. The NPS sensing structures are evenly distributed over the surface.



## Nanoplasmonic sensors

Insplorion's nanoplasmonic sensors are produced in a state of the art clean room facility, which guarantees stability, reproducibility and well controlled surface chemistry. The nanostructures cover the whole sensor area, forming a quasi random pattern on the surface. The structures are uniform in size and shape. On Insplorion's standard sensors, the nanostructures are made of gold but other metals are available upon request. The nanostructures can be ordered uncoated or covered by a thin top coating.

	Top coating	Nanostructure type	Nanostructure material	Temperature range	Plasmon peak position
Standard	SiO <sub>2</sub> , Al <sub>2</sub> O <sub>3</sub> , TiO <sub>2</sub>	Discs, embedded discs (topographically flat surface)	Au	RT-150°C	700-780 nm**
Upon request	E.g. AZO, IOH, iron oxides, porous Al <sub>2</sub> O <sub>3</sub>	Spherical faceted particles*, rings, cones, holes, wells, caves*	Ag, Al, Cu, Pd, Pt	RT-600°C*	500-1100 nm

\* not available for Acoulyte sensors  
\*\* in water



**Insplorion™**

Sahlgrenska Science Park  
Medicinaregatan 8A  
413 90 Göteborg  
Sweden  
[www.insplorion.com](http://www.insplorion.com)

## Standard sensors

### Discs

Disc diameter: 100 nm\*  
Disc height: 20 nm\*\*



## Structures made upon request

### Embedded discs (topographically flat surface)

Disc diameter: 100 nm  
Disc height: 20 nm  
Thickness of silicon dioxide layer  
above the gold nanodiscs: 5-10 nm



### Caves

Hole diameter: 60-160 nm  
Depth: 30-120 nm




### Cones/ truncated cones

Base diameter: 60-200 nm  
Cone height:  $\leq 0.8 \times$  diameter  
Sidewall angle:  $30^\circ$  off normal



### Rings

Outer diameter:  $\sim 50$  nm  
Inner diameter:  $\sim 30$  nm  
Ring height:  $\sim 20$  nm




### Spherical faceted particles

Diameter: 60-200 nm




### Holes

Hole diameter: 60-200 nm  
Hole depth:  $\leq$  hole radius



### Wells

Diameter: 70-160 nm  
Depth: 50-200 nm



\* Diameters 60-200 nm are available upon request.

\*\* Heights 10-30 nm are available upon request.

**Insplorion™**


Sahlgrenska Science Park  
Medicinaregatan 8A  
413 90 Göteborg  
Sweden  
[www.insplorion.com](http://www.insplorion.com)

## Examples

### Influence of surface chemistry on membrane architecture

Insplorion's standard sensors with different top coatings (no coating,  $\text{TiO}_2$  and  $\text{SiO}_2$  respectively) were used to study the influence of surface chemistry on lipid vesicle adsorption and bilayer formation. Vesicle adsorption was observed on all surface, while vesicle rupture and lipid bilayer formation was observed only on  $\text{SiO}_2$ .


See application note #14 for further information.



### Biocorona formation on silica nanoparticles

Spherical faceted nanoparticles have been used for in situ characterisation of molecule-nanoparticle interactions. The nanostructures were coated by  $\text{SiO}_2$  to mimic silica nanoparticles in solution. In this study the role of flat versus curved surfaces on the faceted nanoparticles in protein corona formation was quantitatively analyzed by using spherical structures with different diameters and, consequently, different flat/curved ratios.

See application note #28 for further information.



### Protein adsorption to lipid membranes with negative curvature

Using nanowell structures it is possible to distinguish the signal arising from molecular binding inside the holes from that of molecular binding to the flat top surface. The nanowells have been used to study how the binding of proteins to lipids is affected if the membrane has a negative curvature.

See application note #22 for further information.

